

Anti-NUDT15 antibody (62-112) (STJA0006306)

STJA0006306

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nucleotide triphosphate diphosphatase NUDT15 (62-112) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55270
Gene Symbol	NUDT15
Uniprot ID	NUDT15_HUMAN
Immunogen	Synthesized peptide derived from human NUD15 AA range: 62-112
Immunogen Region	62-112
Specificity	This antibody detects endogenous levels of human NUD15
Immunogen Sequence	