

Anti-PSMB7 antibody (191-241) (STJA0006321)

STJA0006321

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Proteasome subunit beta type-7 (191-241) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5695
Gene Symbol	PSMB7
Uniprot ID	PSB7_HUMAN
Immunogen	Synthesized peptide derived from human PSB7 AA range: 191-241
Immunogen Region	191-241
Specificity	This antibody detects endogenous levels of human PSB7
Immunogen Sequence	