

Anti-OLFM1 antibody (414-464) (STJA0006322)

STJA0006322

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Noelin (414-464) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10439
Gene Symbol	OLFM1
Uniprot ID	NOE1_HUMAN
Immunogen	Synthesized peptide derived from human NOE1 AA range: 414-464
Immunogen Region	414-464
Specificity	This antibody detects endogenous levels of human NOE1
Immunogen Sequence	