

Anti-LHX6 antibody (257-307) (STJA0006324)

STJA0006324

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LIM/homeobox protein Lhx6 (257-307) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26468
Gene Symbol	LHX6
Uniprot ID	LHX6_HUMAN
Immunogen	Synthesized peptide derived from human LHX6 AA range: 257-307
Immunogen Region	257-307
Specificity	This antibody detects endogenous levels of human LHX6
Immunogen Sequence	