

Anti-MBD3L1 antibody (136-186) (STJA0006328)

STJA0006328

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Methyl-CpG-binding domain protein 3-like 1 (136-186) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85509
Gene Symbol	MBD3L1
Uniprot ID	MB3L1_HUMAN
Immunogen	Synthesized peptide derived from human MB3L1 AA range: 136-186
Immunogen Region	136-186
Specificity	This antibody detects endogenous levels of human MB3L1
Immunogen Sequence	