

Anti-ZC3H12D antibody (87-137) (STJA0006332)

STJA0006332

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Probable ribonuclease ZC3H12D (87-137) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	340152
Gene Symbol	ZC3H12D
Uniprot ID	ZC12D_HUMAN
Immunogen	Synthesized peptide derived from human ZC12D AA range: 87-137
Immunogen Region	87-137
Specificity	This antibody detects endogenous levels of human ZC12D
Immunogen Sequence	