

Anti-SERHL2 antibody (259-309) (STJA0006336)

STJA0006336

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine hydrolase-like protein 2 (259-309) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	253190
Gene Symbol	SERHL2
Uniprot ID	SEHL2_HUMAN
Immunogen	Synthesized peptide derived from human SEHL2 AA range: 259-309
Immunogen Region	259-309
Specificity	This antibody detects endogenous levels of human SEHL2
Immunogen Sequence	