

Anti-LAD1 antibody (404-454) (STJA0006338)

STJA0006338

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ladinin-1 (404-454) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3898
Gene Symbol	LAD1
Uniprot ID	LAD1_HUMAN
Immunogen	Synthesized peptide derived from human LAD1 AA range: 404-454
Immunogen Region	404-454
Specificity	This antibody detects endogenous levels of human LAD1
Immunogen Sequence	