

Anti-CHRFAM7A antibody (124-174) (STJA0006350)

STJA0006350

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CHRNA7-FAM7A fusion protein (124-174) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89832
Gene Symbol	CHRFAM7A
Uniprot ID	CRFM7_HUMAN
Immunogen	Synthesized peptide derived from human CRFM7 AA range: 124-174
Immunogen Region	124-174
Specificity	This antibody detects endogenous levels of human CRFM7
Immunogen Sequence	