

Anti-VIRMA antibody (1285-1335) (STJA0006367)

STJA0006367

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein virilizer homolog (1285-1335) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25962
Gene Symbol	VIRMA
Uniprot ID	VIR_HUMAN
Immunogen	Synthesized peptide derived from human VIR AA range: 1285-1335
Immunogen Region	1285-1335
Specificity	This antibody detects endogenous levels of human VIR
Immunogen Sequence	