

Anti-NABP2 antibody (70-120) (STJA0006368)

STJA0006368

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SOSS complex subunit B1 (70-120) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79035
Gene Symbol	NABP2
Uniprot ID	SOSB1_HUMAN
Immunogen	Synthesized peptide derived from human SOSB1 AA range: 70-120
Immunogen Region	70-120
Specificity	This antibody detects endogenous levels of human SOSB1
Immunogen Sequence	