

Anti-ZNF444 antibody (189-239) (STJA0006381)

STJA0006381

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Zinc finger protein 444 (189-239) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55311
Gene Symbol	ZNF444
Uniprot ID	ZNF444_HUMAN
Immunogen	Synthesized peptide derived from human ZN444 AA range: 189-239
Immunogen Region	189-239
Specificity	This antibody detects endogenous levels of human ZN444
Immunogen Sequence	