

Anti-GEMIN4 antibody (754-804) (STJA0006395)

STJA0006395

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Gem-associated protein 4 (754-804) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50628
Gene Symbol	GEMIN4
Uniprot ID	GEMIN4_HUMAN
Immunogen	Synthesized peptide derived from human GEMIN4 AA range: 754-804
Immunogen Region	754-804
Specificity	This antibody detects endogenous levels of human GEMIN4
Immunogen Sequence	