

Anti-RASGRP2 antibody (492-542) (STJA0006437)

STJA0006437

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RAS guanyl-releasing protein 2 (492-542) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10235
Gene Symbol	RASGRP2
Uniprot ID	GRP2_HUMAN
Immunogen	Synthesized peptide derived from human GRP2 AA range: 492-542
Immunogen Region	492-542
Specificity	This antibody detects endogenous levels of human GRP2
Immunogen Sequence	