

Anti-NUMBL antibody (290-340) (STJA0006441)

STJA0006441

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Numb-like protein (290-340) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9253
Gene Symbol	NUMBL
Uniprot ID	NUMBL_HUMAN
Immunogen	Synthesized peptide derived from human NUMBL AA range: 290-340
Immunogen Region	290-340
Specificity	This antibody detects endogenous levels of human NUMBL
Immunogen Sequence	