

Anti-GTPBP6 antibody (410-460) (STJA0006477)

STJA0006477

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Putative GTP-binding protein 6 (410-460) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8225
Gene Symbol	GTPBP6
Uniprot ID	GTPB6_HUMAN
Immunogen	Synthesized peptide derived from human GTPB6 AA range: 410-460
Immunogen Region	410-460
Specificity	This antibody detects endogenous levels of human GTPB6
Immunogen Sequence	