

Anti-PIK3R4 antibody (170-220) (STJA0006495)

STJA0006495

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phosphoinositide 3-kinase regulatory subunit 4 (170-220) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	30849
Gene Symbol	PIK3R4
Uniprot ID	PI3R4_HUMAN
Immunogen	Synthesized peptide derived from human PI3R4 AA range: 170-220
Immunogen Region	170-220
Specificity	This antibody detects endogenous levels of human PI3R4
Immunogen Sequence	