

Anti-FLYWCH2 antibody (44-94) (STJA0006508)

STJA0006508

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FLYWCH family member 2 (44-94) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114984
Gene Symbol	FLYWCH2
Uniprot ID	FWCH2_HUMAN
Immunogen	Synthesized peptide derived from human FWCH2 AA range: 44-94
Immunogen Region	44-94
Specificity	This antibody detects endogenous levels of human FWCH2
Immunogen Sequence	