

Anti-EEF1E1 antibody (70-120) (STJA0006519)

STJA0006519

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Eukaryotic translation elongation factor 1 epsilon-1 (70-120) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9521
Gene Symbol	EEF1E1
Uniprot ID	MCA3_HUMAN
Immunogen	Synthesized peptide derived from human MCA3 AA range: 70-120
Immunogen Region	70-120
Specificity	This antibody detects endogenous levels of human MCA3
Immunogen Sequence	