

Anti-SPINK1 antibody (26-75) (STJA0006545)

STJA0006545

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine protease inhibitor Kazal-type 1 (26-75) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6690
Gene Symbol	SPINK1
Uniprot ID	ISK1_HUMAN
Immunogen	Synthesized peptide derived from human ISK1 AA range: 26-75
Immunogen Region	26-75
Specificity	This antibody detects endogenous levels of human ISK1
Immunogen Sequence	