

Anti-TPM3 antibody (175-225) (STJA0006598)

STJA0006598

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tropomyosin alpha-3 chain (175-225) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7170
Gene Symbol	TPM3
Uniprot ID	TPM3_HUMAN
Immunogen	Synthesized peptide derived from human TPM3 AA range: 175-225
Immunogen Region	175-225
Specificity	This antibody detects endogenous levels of human TPM3
Immunogen Sequence	