

Anti-HOOK2 antibody (224-274) (STJA0006607)

STJA0006607

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Hook homolog 2 (224-274) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29911
Gene Symbol	HOOK2
Uniprot ID	HOOK2_HUMAN
Immunogen	Synthesized peptide derived from human HOOK2 AA range: 224-274
Immunogen Region	224-274
Specificity	This antibody detects endogenous levels of human HOOK2
Immunogen Sequence	