

**Anti-THOC6 antibody (79-129) (STJA0006610)**

STJA0006610

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-THO complex subunit 6 homolog (79-129) is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB  |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Human/Mouse/Rat   |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 1 mg/mL   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:1000-2000  |
| <b>Formulation</b>         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

**TARGET INFORMATION**

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | <a href="#">79228</a>   |
| <b>Gene Symbol</b>        | <a href="#">THOC6</a>   |
| <b>Uniprot ID</b>         | <a href="#">THOC6_HUMAN</a>                                   |
| <b>Immunogen</b>          | Synthesized peptide derived from human THOC6 AA range: 79-129 |
| <b>Immunogen Region</b>   | 79-129  |
| <b>Specificity</b>        | This antibody detects endogenous levels of human THOC6        |
| <b>Immunogen Sequence</b> |   |