

Anti-VBP1 antibody (71-121) (STJA0006625)

STJA0006625

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Prefoldin subunit 3 (71-121) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7411
Gene Symbol	VBP1
Uniprot ID	PFD3_HUMAN
Immunogen	Synthesized peptide derived from human PFD3 AA range: 71-121
Immunogen Region	71-121
Specificity	This antibody detects endogenous levels of human PFD3
Immunogen Sequence	