

Anti-NEU1 antibody (139-189) (STJA0006640)

STJA0006640

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Sialidase-1 (139-189) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4758
Gene Symbol	NEU1
Uniprot ID	NEUR1_HUMAN
Immunogen	Synthesized peptide derived from human NEUR1 AA range: 139-189
Immunogen Region	139-189
Specificity	This antibody detects endogenous levels of human NEUR1
Immunogen Sequence	