

Anti-PDE6C antibody (474-524) (STJA0006643)

STJA0006643

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cone cGMP-specific 3 NA, 5 cyclic phosphodiesterase subunit alpha (474-524) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5146
Gene Symbol	PDE6C
Uniprot ID	PDE6C_HUMAN
Immunogen	Synthesized peptide derived from human PDE6C AA range: 474-524
Immunogen Region	474-524
Specificity	This antibody detects endogenous levels of human PDE6C
Immunogen Sequence	