

Anti-ERN1 antibody (282-433) (STJA0006696)

STJA0006696

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine/threonine-protein kinase/endoribonuclease IRE1 (282-433) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2081
Gene Symbol	ERN1
Uniprot ID	ERN1_HUMAN
Immunogen	Purified recombinant fragment of human IRE1 Alpha (aa282-433) expressed in E. Coli.
Immunogen Region	282-433
Specificity	IRE1 Alpha Antibody detects endogenous levels of IRE1 Alpha protein.
Immunogen Sequence	