

Anti-GAPDH-Drosophilidae antibody (Full Length) (STJA0006712)

STJA0006712

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GAPDH-Drosophilidae (Full Length) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	D.melanogaster

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified recombinant Full length GAPDH (Drosophilidae)
Immunogen Region	Full Length
Specificity	This antibody detects endogenous levels of GAPDH (Drosophilidae)
Immunogen Sequence	



Western Blot analysis of Drosophilidae cells using GAPDH rabbit pAb (Drosophilidae). Secondary antibody (STJ99566) was diluted at 1:20000.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081