

Anti-STAU1 antibody (STJA0006849)

STJA0006849

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Double-stranded RNA-binding protein Staufen homolog 1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6780
Gene Symbol	STAU1
Uniprot ID	STAU1_HUMAN
Immunogen	Synthesized peptide derived from human STAU1
Immunogen Region	
Specificity	This antibody detects endogenous levels of STAU1 at Human, Mouse
Immunogen Sequence	