

Anti-VTCN1 antibody (STJA0007171)

STJA0007171

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-V-set domain-containing T-cell activation inhibitor 1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79679
Gene Symbol	VTCN1
Uniprot ID	VTCN1_HUMAN
Immunogen	Synthesized peptide derived from human VTCN1
Immunogen Region	
Specificity	This antibody detects endogenous levels of VTCN1 at Human, Mouse, Rat
Immunogen Sequence	