

Anti-ECH1 antibody (STJA0007320)

STJA0007320

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Delta3, 5-Delta2, 4-dienoyl-CoA isomerase, mitochondrial is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1891
Gene Symbol	ECH1
Uniprot ID	ECH1_HUMAN
Immunogen	Synthesized peptide derived from human ECH1
Immunogen Region	
Specificity	This antibody detects endogenous levels of ECH1 at Human, Mouse, Rat
Immunogen Sequence	