

Anti-SMPDL3B antibody (STJA0007489)

STJA0007489

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acid sphingomyelinase-like phosphodiesterase 3b is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27293
Gene Symbol	SMPDL3B
Uniprot ID	ASM3B_HUMAN
Immunogen	Synthesized peptide derived from human SMPDL3B
Immunogen Region	
Specificity	This antibody detects endogenous levels of SMPDL3B at Human, Mouse
Immunogen Sequence	