

Anti-ST3GAL4 antibody (STJA0007493)

STJA0007493

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CMP-N-acetylneuraminate-beta-galactosamide-alpha-2, 3-sialyltransferase 4 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6484
Gene Symbol	ST3GAL4
Uniprot ID	SIA4C_HUMAN
Immunogen	Synthesized peptide derived from human ST3GAL4
Immunogen Region	
Specificity	This antibody detects endogenous levels of ST3GAL4 at Human, Mouse, Rat
Immunogen Sequence	