

Anti-CTR9 antibody (1150-1173) (STJA0007522)

STJA0007522

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RNA polymerase-associated protein CTR9 homolog (1150-1173) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9646
Gene Symbol	CTR9
Uniprot ID	CTR9_HUMAN
Immunogen	Synthesized peptide derived from human CTR9.AA range:1150-1173
Immunogen Region	1150-1173
Specificity	This antibody detects endogenous levels of CTR9 at Human, Mouse
Immunogen Sequence	