

Anti-TIMM8A antibody (STJA0007534)

STJA0007534

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Mitochondrial import inner membrane translocase subunit Tim8 A is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1678
Gene Symbol	TIMM8A
Uniprot ID	TIM8A_HUMAN
Immunogen	Synthesized peptide derived from human TIMM8A
Immunogen Region	
Specificity	This antibody detects endogenous levels of TIMM8A at Human, Mouse, Rat
Immunogen Sequence	