

Anti-BLOC1S5 antibody (STJA0007655)

STJA0007655

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Biogenesis of lysosome-related organelles complex 1 subunit 5 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	63915
Gene Symbol	BLOC1S5
Uniprot ID	BL1S5_HUMAN
Immunogen	Synthesized peptide derived from human MUTED/BLOC1S5
Immunogen Region	
Specificity	This antibody detects endogenous levels of MUTED/BLOC1S5 at Human
Immunogen Sequence	