

Anti-ABHD17A antibody (STJA0007854)

STJA0007854

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Alpha/beta hydrolase domain-containing protein 17A is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81926
Gene Symbol	ABHD17A
Uniprot ID	AB17A_HUMAN
Immunogen	Synthesized peptide derived from human ABHD17A
Immunogen Region	
Specificity	This antibody detects endogenous levels of ABHD17A at Human
Immunogen Sequence	