

Anti-TMEM214 antibody (STJA0008082)

STJA0008082

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Transmembrane protein 214 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54867
Gene Symbol	TMEM214
Uniprot ID	TM214_HUMAN
Immunogen	Synthesized peptide derived from human TMEM214
Immunogen Region	
Specificity	This antibody detects endogenous levels of TMEM214 at Human, Mouse, Rat
Immunogen Sequence	