

Anti-APMAP antibody (STJA0008170)

STJA0008170

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Adipocyte plasma membrane-associated protein is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range | WB 1:500-2000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 57136 |
| Gene Symbol | APMAP |
| Uniprot ID | APMAP_HUMAN |
| Immunogen | Synthesized peptide derived from human C20orf3/APMAP |
| Immunogen Region | |
| Specificity | This antibody detects endogenous levels of C20orf3/APMAP at Human, Mouse, Rat |
| Immunogen Sequence | |