

Anti-VPS4A antibody (400-437) (STJA0008215)

STJA0008215

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Vacuolar protein sorting-associated protein 4A (400-437) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27183
Gene Symbol	VPS4A
Uniprot ID	VPS4A_HUMAN
Immunogen	Synthesized peptide derived from human VPS4A. AA range:400-437
Immunogen Region	400-437
Specificity	This antibody detects endogenous levels of VPS4A at Human, Mouse, Rat
Immunogen Sequence	