

Anti-OSGEPL1 antibody (STJA0008265)

STJA0008265

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-tRNA N6-adenosine threonylcarbamoyltransferase, mitochondrialA synthaseA37 threonylcarbamoyladenine biosynthesis protein OSGEPL1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64172
Gene Symbol	OSGEPL1
Uniprot ID	OSGL1_HUMAN
Immunogen	Synthesized peptide derived from human OSGEPL1
Region	
Specificity	This antibody detects endogenous levels of OSGEPL1 at Human, Mouse, Rat
Immunogen Sequence	