

Anti-PIP4K2A antibody (STJA0008345)

STJA0008345

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phosphatidylinositol 5-phosphate 4-kinase type-2 alphaP 4-kinase type II alphaP-5-kinase B isoformP-5-kinase C isoformP-4-kinase isoform 2-alpha is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5305
Gene Symbol	PIP4K2A
Uniprot ID	PI42A_HUMAN
Immunogen	Synthesized peptide derived from human PIP4K2A
Immunogen Region	
Specificity	This antibody detects endogenous levels of PIP4K2A at Human, Mouse, Rat
Immunogen Sequence	