

Anti-IVNS1ABP antibody (STJA0008407)

STJA0008407

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Influenza virus NS1A-binding protein is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10625
Gene Symbol	IVNS1ABP
Uniprot ID	NS1BP_HUMAN
Immunogen	Synthesized peptide derived from human IVNS1ABP
Immunogen Region	
Specificity	This antibody detects endogenous levels of IVNS1ABP at Human, Mouse
Immunogen Sequence	