

Anti-DNAAF6 antibody (STJA0008525)

STJA0008525

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dynein axonemal assembly factor 6 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	139212
Gene Symbol	DNAAF6
Uniprot ID	DAAF6_HUMAN
Immunogen	Synthesized peptide derived from human PIH1D3
Immunogen Region	
Specificity	This antibody detects endogenous levels of PIH1D3 at Human
Immunogen Sequence	