

Anti-ACYP1 antibody (STJA0008568)

STJA0008568

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acylphosphatase-1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	97
Gene Symbol	ACYP1
Uniprot ID	ACYP1_HUMAN
Immunogen	Synthesized peptide derived from human ACYP1
Immunogen Region	
Specificity	This antibody detects endogenous levels of ACYP1 at Human, Mouse
Immunogen Sequence	