

Anti-ZNF382 antibody (STJA0008647)

STJA0008647

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Zinc finger protein 382 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84911
Gene Symbol	ZNF382
Uniprot ID	ZN382_HUMAN
Immunogen	Synthesized peptide derived from human ZNF382
Immunogen Region	
Specificity	This antibody detects endogenous levels of ZNF382 at Human, Mouse, Rat
Immunogen Sequence	