

Anti-BHMT2 antibody (STJA0008671)

STJA0008671

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-S-methylmethionine--homocysteine S-methyltransferase BHMT2 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23743
Gene Symbol	BHMT2
Uniprot ID	BHMT2_HUMAN
Immunogen	Synthesized peptide derived from human BHMT2
Immunogen Region	
Specificity	This antibody detects endogenous levels of BHMT2 at Human, Mouse, Rat
Immunogen Sequence	