

Anti-CAMK1D antibody (STJA0008681)

STJA0008681

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Calcium/calmodulin-dependent protein kinase type 1D is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57118
Gene Symbol	CAMK1D
Uniprot ID	KCC1D_HUMAN
Immunogen	Synthesized peptide derived from human CAMK1D
Immunogen Region	
Specificity	This antibody detects endogenous levels of CAMK1D at Human, Mouse
Immunogen Sequence	