

Anti-ANKRD13A antibody (STJA0008742)

STJA0008742

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ankyrin repeat domain-containing protein 13A is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	88455
Gene Symbol	ANKRD13A
Uniprot ID	AN13A_HUMAN
Immunogen	Synthesized peptide derived from human ANKRD13A
Immunogen Region	
Specificity	This antibody detects endogenous levels of ANKRD13A at Human, Mouse
Immunogen Sequence	