

Anti-SYNE4 antibody (STJA0008773)

STJA0008773

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nesprin-4 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	163183
Gene Symbol	SYNE4
Uniprot ID	SYNE4_HUMAN
Immunogen	Synthesized peptide derived from human Nesprin 4
Immunogen Region	
Specificity	This antibody detects endogenous levels of Nesprin 4 at Human, Mouse, Rat
Immunogen Sequence	